


LC/MS Application Note

On-line automated protein precipitation preparation followed by LC-MS/MS analysis and LDTD-MS/MS cross validation



LIFE SCIENCE



www.palsystem.com

On-line automated protein precipitation preparation followed by LC-MS/MS analysis and LDTD-MS/MS cross validation

Pascal Belisle, Sylvain Letarte, Serge Auger and Pierre Picard
Phytronix Technologies, Quebec, Canada

Introduction

To minimize contamination or inversion during the preanalysis sample preparation, a new on-line sample preparation system is evaluated. The Robotic Tool Change (RTC) is an on-line extraction and injection system. Complex automated sample preparation steps such as dilution, derivatization, mixing, liquid/liquid extraction, and protein precipitation can be performed prior to injection. Using the RTC system, an on-line protein precipitation of plasma samples followed by LC-MS/MS analysis was performed. During the same process, samples were also deposited in LazWell plates for Laser Diode Thermal Desorption analysis (LDTD-MS/MS). The same on-line sample preparation steps were used to do cross validation of LC-MS/MS versus High throughput LDTD-MS/MS analysis. The on-line extraction of Dextrophan / Dextromethorphan plasma samples was used for the cross validation method test.

RTC sample preparation method

Plasma samples (100 μ L) were put in injection vials, capped then transferred to the RTC system. The on-line protein precipitation was performed as follows: automated addition of 800 μ L acetonitrile solution containing internal standard in sample vial using 1 mL syringe and vial was transferred to a vortex mixer. During the phase separation, RTC syringe tool automatically changes from a 1 mL to a 10 μ L and acetonitrile phase is injected in LC-MS/MS system followed by LazWell plate spotting. LC-MS/MS and LDTD-MS/MS analysis were performed on a Sciex QTrap 5500.

LDTD-APCI-MS/MS method

For the LDTD-MS/MS analysis, following the on-line extraction, 5 μ L of protein precipitation sample was added in LazWell plate and evaporated to dryness. Samples were desorbed using the following laser pattern: increase laser power to 25% in 6 seconds, maintain 2 seconds at 25% laser power and decrease laser power to 0% in 0.1 second. An air flow of 3 L/min is used as a carrier gas. Positive APCI ionization mode with 3 μ A discharge current was used with 258 \rightarrow 199 and 272 \rightarrow 215 as transitions for Dextrophan and Dextromethorphan respectively.

LC-MS/MS method

After the on-line extraction, 5 μ L of protein precipitation sample was injected on a SB-C18 column (4.6 X 200, 5 μ m) for LC-MS/MS analysis. A flow rate of 1 mL/min (isocratic mode) using Acetonitrile:Water (50:50) with 0.1% FA as mobile phase was used. Drugs were detected in positive ESI mode using 258 \rightarrow 157 and 272 \rightarrow 215 as transition for Dextrophan (DEX) and Dextromethorphan (DMX) respectively. The retention time of Dextrophan and Dextromethorphan were 3.2 and 4.5 minutes respectively.



Fig. 1: RTC system and a LC coupled to a MS/MS Sciex QTRAP[®] 5500 System in MRM mode

Results and Discussion

LC-MS/MS

Linearity of 0.99638 (DEX) and 0.99842 (DMX) with a calibration range of 0.5 to 100 ng/mL is obtained (Figures 2 & 3). Intra-run accuracy is between 93 and 115 % and precision is below 11 % (Table 1).




Fig. 2: DMX calibration curve with LC-MS/MS



Fig. 3: DEX calibration curve with LC-MS/MS

Compound	Conc (ng/mL)	n	Mean (ng/mL)	%RSD	%NOM
DEX	0.5	4	0.57	6.22	114.76
	2.5		2.33	6.43	93.14
	10		9.72	5.21	97.20
	50		50.33	6.52	100.67
	100		100.80	9.45	100.80
DMX	0.5	4	0.53	10.36	106.39
	2.5		2.39	8.41	95.78
	10		9.82	7.83	98.17
	50		50.43	4.53	100.86
	100		99.48	6.50	99.48

Tab. 1: DEX and DMX intra run accuracy with LC-MS/MS

LDTD-MS/MS

Linearity of 0.99785 (DEX) and 0.99791 (DMX) with a calibration range of 0.5 to 100 ng/mL is obtained (Figures 4 & 5). Intra-run accuracy is between 80 and 107 % and precision is below 10% (Table 2).




Fig. 4: DMX calibration curve with LDTD-MS/MS




Fig. 5: DEX calibration curve with LDTD-MS/MS

Cross validation

Cross validation of unknown samples was performed (Figures 6 & 7). A correlation of 0.9987 and 0.9981 is obtained between LC-MS/MS and LDTD-MS/MS results for Dextrorphan and Dextromethorphan, respectively.

Compound	Conc (ng/mL)	n	Mean (ng/mL)	%RSD	%NOM
DEX	0.5	4	0.43	8.34	86.13
	2.5		2.48	1.13	99.24
	10		10.61	1.85	106.06
	50		52.40	1.39	104.81
	100		94.70	3.87	94.70
DMX	0.5	4	0.41	8.94	81.08
	2.5		2.54	3.95	101.47
	10		10.69	1.10	106.87
	50		52.44	1.71	104.88
	100		94.48	1.67	94.48

Tab. 2: DEX and DMX intra run accuracy with LDTD-MS/MS




Fig. 6: Cross-validation between LC and LDTD (Dextrorphan)




Fig. 7: Cross-validation between LC and LDTD (Dextromethorphan)

Conclusions

The Robotic Tool Change (RTC) effectively performs on-line automated protein precipitation for the analysis of Dextromethorphan and Dextrorphan. The extracts were analysed using 2 techniques: LC-MS/MS and LDTD-MS/MS. The RTC is an excellent tool for a cross-validation between two different techniques.

For more information visit www.phytronics.com

Legal Statements

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described in this document at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall CTC Analytics AG be held liable for any incidental damage or damages arising as a consequence of or from the use of this document.

© 2016 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Imprint

Date of print: 04.2016

CTC Analytics AG
Industriestrasse 20
CH-4222 Zwingen
Switzerland
T +41 61 765 81 00
Contact: info@ctc.ch