


## Pyrolysis of Rubber with Antioxidant 6-PPD

Rubber products are complex polymeric materials which generally include copolymers, fillers, antioxidants, plasticisers and other additives. Automobile tire rubbers are good examples, since they usually contain significant amounts of carbon black, hydrocarbon oil and antioxidants in a matrix that may include styrene, butadiene and isoprene. Nevertheless, a single pyrolysis run can provide significant information about such a complex material.


Figure 1 shows the pyrolysis at 700°C of a sample of tire rubber. The matrix is an isoprene/butadiene copolymer, and early in the pyrogram the monomers and dimers from the copolymer may be seen. Later in the chromatogram is a peak for the antioxidant 6-PPD [N-(1,3-dimethyl butyl)-N'-phenyl-p-phenylenediamine]. 6-PPD has a molecular weight of 268, and produces the mass spectrum seen in Figure 2.

Other antioxidants, especially hindered phenolics, are also frequently seen in the pyrolysis of rubber samples. Some of these, like BHT, are present as the parent molecule. Other antioxidant molecules are too large to permit GC analysis intact, and are seen as fragments, just as for the polymer itself.

At the end of the pyrogram is a region of unresolved peaks (marked 6 in Figure 1) which results from the oil added to the rubber in making the tire. Consequently, from just a single pyrolysis run, the monomer composition of the rubber, the antioxidant used and the oil added may be determined.



**Figure 1.** Pyrolysis of tire rubber at 700°C. Peak 1=butadiene, 2=isoprene, 3=butadiene dimer, 4=isoprene dimer, 5= 6-PPD, 6=Hydro-carbon oil.



**Figure 2.** Mass spectrum of peak at 21.5 min.

## Equipment

The rubber sample was analyzed using a CDS Pyroprobe Model 2000, interfaced to an Agilent 6890 gas chromatograph which was equipped with an Model 5973 MSD as the detector. A 100µg piece of the rubber was placed into a quartz tube, which was inserted into the coil filament of the Pyroprobe.

### Pyroprobe 2000 Conditions

Interface: Model 1500 valved  
Interface temp: 300°C isothermal 700°  
Temperature: C for 10 seconds 10°  
Rate: C/ms  
Probe: Platinum coil

### GC Conditions

Carrier: Helium  
Split: 50:1  
Column: HP-5 (30m X 0.25)  
Detector: MSD  
GC Program:  
Initial: 40°C for 2 minutes  
Ramp: 10°C/min.  
Final: 300°C

FOR MORE INFORMATION  
CONCERNING THIS APPLICATION,  
WE RECOMMEND THE  
FOLLOWING READING:

Temperature as a Sample Preparation  
Tool in the Analysis of Materials by  
GC-MS, T. P. Wampler, LC.GC, 17  
(1999) S14.

© courtesy of CDS Analytical Inc.

## INTERNATIONAL

JSB Group  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 251 47 53 F  
+31 (0) 40 251 47 58

[INFO@GO-JSB.COM](mailto:INFO@GO-JSB.COM)  
[WWW.GO-JSB.COM](http://WWW.GO-JSB.COM)

## SALES AND SERVICE

### NETHERLANDS

Apolloweg 2B  
8239 DA Lelystad  
T +31 (0) 32 087 00 18  
F +31 (0) 32 087 00 19

### GERMANY, AUSTRIA, SWITZERLAND

Max-Planck-Strasse 4  
D-47475 Kamp-Lintfort  
T +49 (0) 28 42 9280 799  
F +49 (0) 28 42 9732 638

### BELGIUM

Grensstraat 7, Box 3  
1831 Diegem  
T +32 (0) 27219211  
F +32 (0) 27207622

### UK & IRELAND

Cedar Court, Grove Park Business Est. White  
Waltham, Maidenhead, Berks, SL6 3LW T +44  
(0) 16 288 220 48  
F +44 (0) 70 394 006 78

