



## Phthalate Calibration and Reproducibility using an IEC Standard Method

Application Note  
Electronics Industry

Author:

Karen Sam

### Abstract

This application note presents RSDs of phthalates for IEC 62321-8 method for using a CDS 6000 Series Pyroprobe Autosampler.

Materials being investigated by thermal sampling techniques such as pyrolysis are frequently polymeric, but may also have volatile and semi-volatile contaminants and additives. When using pyrolysis- gas chromatography, these non-polymeric constituents produce some of the most significant peaks in the chromatogram. It has become common to perform two or more analyses on the sample at increasing temperatures, to remove the additives before pyrolysis. In some instances, however, analysis of the polymer structure takes a back seat to the nature and amount of additives present.

The International Electrotechnical Commission (IEC) published a new standard for hazardous substances in electronic equipment for determining phthalates in polymeric materials. IEC 62321-8 defines approaches to determine phthalates di-isobutyl phthalate (DIBP), di-n-butyl phthalate (DBP), benzylbutyl phthalate (BBP), bis-2-ethyl hexyl phthalate (DEHP), di-n-octyl phthalate (DNOP), di-isononyl phthalate (DINP) and di-iso-decyl phthalate (DIDP) in electronics, by GC-MS, and TD-GC-MS. The method involves two separate heating ramps for one GC run. This can be performed using a CDS 6000 Series Pyroprobe Autosampler. TIC and extracted ion chromatograms in Figure 1 match the chromatograms in Annex C.2 of the International Standard.



The diagnostic worth of the results depends on the reproducibility, and reproducibility for thermal desorption depends greatly on temperature precision, along with sample related issues like homogeneity and sample preparation. Like all analytical instrumentation, the Pyroprobe is designed to perform with optimum precision. Pyroprobe filaments are calibrated using optical pyrometry. Using this technique, a series of 20 firings at 1100°C produced an average measured temperature of 1100.15°C with a relative standard deviation of only 0.04%. This reproducibility is not only seen in peak area ratios with pyrolysis of polymers, but also in thermal desorption of phthalate standards. Here, five microliters of a 100ng/mL of a phthalate solution in hexane was added to DISC tubes for replicate analysis in accordance with the IEC method for Phthalates. Figure 2 shows replicate TICs for two of the phthalates, DIBP and DBP. Eight replicates of the standard presented area RSDs for most of the phthalates around or under 3% (Table 1).

### Experimental Parameters

The sample was pyrolyzed in a DISC tube, using a CDS Pyroprobe 6200 with Autosampler.

#### Method 1:

Pyroprobe :

Initial: 200°C  
Ramp: 20°C/minute  
Final: 300°C

Interface:

Rest: 300°C  
Initial: 300°C  
Transfer Line: 300°C  
Valve Oven: 300°C

GC Signal:

GC ready: ON  
GC start: ON

#### Method 2:

Pyroprobe :

Initial: 300°C  
Ramp: 5°C/minute  
Final: 340°C hold 1 min

Interface:

Rest: 300°C  
Initial: 300°C  
Transfer Line: 300°C  
Valve Oven: 300°C

GC Signal:

GC ready: OFF  
GC start: OFF

These two methods were run in sequence during one GC run.

GC/MS

Column: 5% phenyl (30m x 0.25mm)  
Carrier: Helium, 50:1 split  
Injector: 320°C  
Oven: 80°C for 13 minutes  
20°C/min to 300°C  
hold 5 minutes  
Ion Source: 230°C  
Mass Range: 50-1000amu




Figure 1: 500ng Phthalate Standard TIC and Extracted Ion Chromatograms.





Figure 2: Extracted Ion 223 Overlap of 4 Replicates of DIBP and DBP to show similarity of peak areas.

Phthalate	Quant Ion	Area RSD
DIBP	223	3.2 %
DBP	223	2.3 %
BBP	206	4.3 %
DEHP	279	2.9 %
DNOP	279	3.2 %
DINP	293	3.0 %
DIDP	307	3.2 %

Table 1: Area RSDs of 7 regulated phthalates.



Calibration and determination of phthalate concentration is based on a one-point calibration, the area count of each phthalate plotted against its amount. This calibration plot is shown in Figure 3.

The latest version of the Pyroprobe from CDS Analytical ensures repeatable, reliable results for thermal desorption of phthalates in accordance with standard methods, like IEC 6321-8 for determination of phthalates in electrotechnical products.

Figure 3: One Point Calibration for Phthalates.