


User Interface of PAL P&T Plugin



New PAL P&T Plugin Integrates the CDS 7000C Purge and Trap (P&T)

P&T reaches trace sensitivity (ppt) of volatile compounds in water.

Until now, P&T often runs on complex and inflexible setups.

The new PAL System P&T Plugin now enables full control:

- Sampling, extraction, and analysis in one sample list
- Transparent data processing
- Handshake via Plugin: Remote I/O Cables are obsolete
- Flexibility with any other GC sample technique (Headspace, SPME, Liquid)
- Automated Calibration, e.g. with EPA MegaMix® Standards
- Fully meets and exceeds all EPA requirements